

NEB #	Buffers	Size	Stockroom #
B7021S	Gel Loading Dye, Blue (6X)	4 ml	B7021S
NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Stockroom #
E2611S	Gibson Assembly® Master Mix	10 reactions	E2611S
E0552S	Q5® Site-Directed Mutagenesis Kit (Without Competent Cells)	10 reactions	E0552S
NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
E3005S	Luna® Universal One-Step RT-qPCR Kit	200 reactions	E3005S
E3006S	Luna® Universal Probe One-Step RT-qPCR Kit	200 reactions	E3006S
M3004S	Luna® Universal Probe qPCR Master Mix	200 reactions	M3004S
M3003S	Luna® Universal qPCR Master Mix	200 reactions	M3003S
M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	M0530S
M0530L	Phusion™ High-Fidelity DNA Polymerase	500 units	M0530L
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units	M0267S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0174S	AvrII	100 units	R0174S
R3136S	BamHI-HF®	10,000 units	R3136S
R0197S	Clal	1,000 units	R0197S
R0176S	DpnI	1,000 units	R0176S
R3101S	EcoRI-HF®	10,000 units	R3101S
R3104S	HindIII-HF®	10,000 units	R3104S
R3142S	KpnI-HF®	4,000 units	R3142S
R3193S	NcoI-HF®	1,000 units	R3193S
R0111S	NdeI	4,000 units	R0111S
R3131S	NheI-HF®	1,000 units	R3131S
R3189S	NotI-HF®	500 units	R3189S
R0547S	PacI	250 units	R0547S

R3140S	PstI-HF [®]	10,000 units	R3140S
R3150S	PvuI-HF [®]	500 units	R3150S
R0167S	RsaI	1,000 units	R0167S
R3156S	SacI-HF [®]	2,000 units	R3156S
R3138S	Sall-HF [®]	2,000 units	R3138S
R0141S	SmaI	2,000 units	R0141S
R3182S	SphI-HF [®]	500 units	R3182S
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0180S	XmaI	500 units	R0180S

NEB Representative: James DeBoe
jdeboe@neb.com
(978) 778-6011

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com