

NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M2200L	Quick Ligation™ Kit	150 reactions	M2200L
M2200S	Quick Ligation™ Kit	30 reactions	M2200S
M0202T	T4 DNA Ligase	20,000 units	M0202T
M0202S	T4 DNA Ligase	20,000 units	M0202S
NEB #	Epigenetics	Size	Stockroom #
R0171S	HpaII	2,000 units	R0171S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3200L	1 kb Plus DNA Ladder	1,000 gel lanes (0.5 ml)	N3200L
N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	N3200S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
M0210S	DNA Polymerase I, Large (Klenow) Fragment	200 units	M0210S
M0203S	T4 DNA Polymerase	150 units	M0203S
M0267L	Taq DNA Polymerase with ThermoPol® Buffer	2,000 units	M0267L
M0267S	Taq DNA Polymerase with ThermoPol® Buffer	400 units	M0267S
M0254S	Vent® DNA Polymerase	200 units	M0254S
NEB #	Restriction Endonucleases	Size	Stockroom #
R0117S	AatII	500 units	R0117S
R0161S	AccI	1,000 units	R0161S
R0114S	ApaI	5,000 units	R0114S
R0152S	AvaI	2,000 units	R0152S
R3136S	BamHI-HF®	10,000 units	R3136S
R0144S	BglII	2,000 units	R0144S
R0199S	BssHII	500 units	R0199S
R0113S	BstXI	1,000 units	R0113S
R0197S	Clal	1,000 units	R0197S
R3101S	EcoRI-HF®	10,000 units	R3101S
R3195S	EcoRV-HF®	4,000 units	R3195S
R0108S	HaeIII	3,000 units	R0108S
R0104S	HindIII	10,000 units	R0104S
R0104T	HindIII	10,000 units	R0104T

R0105S	HpaI	500 units	R0105S
R3142S	KpnI-HF <sup>®</sup>	4,000 units	R3142S
R0534S	MscI	250 units	R0534S
R3193S	NcoI-HF <sup>®</sup>	1,000 units	R3193S
R0111S	NdeI	4,000 units	R0111S
R3131S	NheI-HF <sup>®</sup>	1,000 units	R3131S
R3189S	NotI-HF <sup>®</sup>	500 units	R3189S
R0127S	NsiI	1,000 units	R0127S
R3140S	PstI-HF <sup>®</sup>	10,000 units	R3140S
R3156S	SacI-HF <sup>®</sup>	2,000 units	R3156S
R0157S	SacII	2,000 units	R0157S
R3138S	SalI-HF <sup>®</sup>	2,000 units	R3138S
R0169S	Sau3AI	200 units	R0169S
R3122S	ScaI-HF <sup>®</sup>	1,000 units	R3122S
R0141S	SmaI	2,000 units	R0141S
R3182S	SphI-HF <sup>®</sup>	500 units	R3182S
R0145S	XbaI	3,000 units	R0145S
R0146S	XhoI	5,000 units	R0146S
R0194S	XmnI	1,000 units	R0194S

**NEB Representative: Mike Mancia**  
mmancia@neb.com  
(978) 380-3941

**NEBnow Support: 1-800-632-5227**  
Email: freezers@neb.com