

Tray	NEB #	Buffers	Size	Quantity Available
13	B7024S	Gel Loading Dye, Purple (6X)	4 ml	1
Tray	NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Quantity Available
8	E2621S	NEBuilder® HiFi DNA Assembly Master Mix	10 reactions	3
Tray	NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Quantity Available
10	M0202S	T4 DNA Ligase	20,000 units	1
10	M0202L	T4 DNA Ligase	100,000 units	1
Tray	NEB #	Glycobiology	Size	Quantity Available
13	P0705S	PNGase F (Glycerol-free)	15,000 units	2
Tray	NEB #	Markers & Ladders	Size	Quantity Available
11	N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	2
11	N3200S	1 kb Plus DNA Ladder	200 gel lanes (0.1 ml)	2
11	N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	1
12	N0552S	Quick-Load® Purple 1 kb DNA Ladder	125 gel lanes (1.25 ml)	2
12	N0550S	Quick-Load® Purple 1 kb Plus DNA Ladder	250 gel lanes (1.25 ml)	1
Tray	NEB #	Nucleic Acid Purification	Size	Quantity Available
SH 1	T2010S	Monarch® Total RNA Miniprep Kit	50 preps	1
SH 4	T1120S	Monarch® Spin DNA Gel Extraction Kit	50 preps	1
SH 1	T2110S	Monarch® Spin RNA Isolation Kit (Mini)	50 preps	1
Tray	NEB #	PCR, Polymerases & Amplification Technologies	Size	Quantity Available
7	N0447S	Deoxynucleotide (dNTP) Solution Mix	8 µmol	1
6	M0530S	Phusion™ High-Fidelity DNA Polymerase	100 units	3
12	M0530L	Phusion™ High-Fidelity DNA Polymerase	500 units	1
6	M0491L	Q5® High-Fidelity DNA Polymerase	500 units	1
5	M0273L	Taq DNA Polymerase with Standard Taq Buffer	2,000 units	1
5	M0273S	Taq DNA Polymerase with Standard Taq Buffer	400 units	1
Tray	NEB #	Protein Tools	Size	Quantity Available
13	P8107S	Proteinase K, Molecular Biology Grade	2 ml	2
Tray	NEB #	Restriction Endonucleases	Size	Quantity Available
2	R0176L	DpnI	5,000 units	1

2	R0176S	Dpnl	1,000 units	out of stock
2	R3101S	EcoRI-HF [®]	10,000 units	1
3	R3104S	HindIII-HF [®]	10,000 units	2
4	R3142S	KpnI-HF [®]	4,000 units	2
3	R3189S	NotI-HF [®]	500 units	1
3	R3156S	SacI-HF [®]	2,000 units	1
4	R0145S	XbaI	3,000 units	2
4	R0146S	XhoI	5,000 units	1

NEB Representative: Allan Waite
awaite@neb.com
(978) 471-0698

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com