

NEB #	Competent Cells	Size	Stockroom #
C2987H	NEB [®] 5-alpha Competent <i>E. coli</i> (High Efficiency)	20 x 0.05 ml	C2987H
C3040H	NEB [®] Stable Competent <i>E. coli</i> (High Efficiency)	20 x 0.05 ml	C3040H
NEB #	DNA Assembly, Cloning and Mutagenesis Kits	Size	Stockroom #
E2611S	Gibson Assembly [®] Master Mix	10 reactions	E2611S
M0554S	KLD Enzyme Mix	25 reactions	M0554S
E5520S	NEBuilder [®] HiFi DNA Assembly Cloning Kit	10 reactions	E5520S
E2621S	NEBuilder [®] HiFi DNA Assembly Master Mix	10 reactions	E2621S
E0554S	Q5 [®] Site-Directed Mutagenesis Kit	10 reactions	E0554S
NEB #	DNA Modifying Enzymes and Cloning Technologies	Size	Stockroom #
M0525S	Quick CIP	1,000 units	M0525S
M2200S	Quick Ligation [™] Kit	30 reactions	M2200S
M0202L	T4 DNA Ligase	100,000 units	M0202L
M0202S	T4 DNA Ligase	20,000 units	M0202S
M0201S	T4 Polynucleotide Kinase	500 units	M0201S
NEB #	Markers & Ladders	Size	Stockroom #
N3232S	1 kb DNA Ladder	200 gel lanes (0.2 ml)	N3232S
N3231S	100 bp DNA Ladder	100 gel lanes (0.1 ml)	N3231S
NEB #	PCR, Polymerases & Amplification Technologies	Size	Stockroom #
N0447S	Deoxynucleotide (dNTP) Solution Mix	8 μmol	N0447S
M0530S	Phusion [™] High-Fidelity DNA Polymerase	100 units	M0530S
M0530L	Phusion [™] High-Fidelity DNA Polymerase	500 units	M0530L
M0494S	Q5 [®] Hot Start High-Fidelity 2X Master Mix	100 reactions	M0494S
M0491L	Q5 [®] High-Fidelity DNA Polymerase	500 units	M0491L
M0491S	Q5 [®] High-Fidelity DNA Polymerase	100 units	M0491S
M0273L	Taq DNA Polymerase with Standard Taq Buffer	2,000 units	M0273L
NEB #	Protein Expression & Purification Technologies	Size	Stockroom #
C2527H	BL21(DE3) Competent <i>E. coli</i>	20 x 0.05 ml	C2527H
NEB #	Protein Tools	Size	Stockroom #
P8107S	Proteinase K, Molecular Biology Grade	2 ml	P8107S
NEB #	Restriction Endonucleases	Size	Stockroom #

R3136S	BamHI-HF [®]	10,000 units	R3136S
R0176S	DpnI	1,000 units	R0176S
R3101S	EcoRI-HF [®]	10,000 units	R3101S
R3193S	NcoI-HF [®]	1,000 units	R3193S
R3131S	NheI-HF [®]	1,000 units	R3131S
R3189L	NotI-HF [®]	2,500 units	R3189L
R3189S	NotI-HF [®]	500 units	R3189S
R0146S	XhoI	5,000 units	R0146S

NEB Representative: James DeBoe
jdeboe@neb.com
(978) 778-6011

NEBnow Support: 1-800-632-5227
Email: freezers@neb.com